

(Demo) Product Specification Cat. No. AG112

Description: Immunoaffinity purification of IgY from egg yolk
 Animal Designation: #106
 Column Matrix: Sepharose – Peptide 14-032 (column 14-069)
 Final Buffer: Tris-Glycine-Saline pH 8

Lot Number :	14-088-B	14-088-D
Immune Status:	Pre-immune	Immune
Number of eggs:	3	7
Volume:	10 ml	12 ml
Protein Concentration:	0.048 mg/ml ¹	0.86 mg/ml ¹
Total Yield:	0.48 milligrams	10.3 milligrams
¹ Protein concentration is determined in a modified dye binding assay using Bovine Gamma Globulin (BGG) as a standard		

