

(Demo) Product Specification Cat. No. AG080

Description: Custom IgY purification from egg yolk
 Animal Designation: #11326
 Immunogen: KLH-Peptide
 Buffer: Phosphate Buffered Saline pH 7.4

Lot Number :	14-054-I	14-054-J
Immune Status:	Pre-immune	Immune
Number of eggs:	3	6
Volume:	20 ml	50 ml
Protein Concentration:	5.63 mg/ml ¹	22.1 mg/ml ¹
Total Yield:	113 milligrams	1105 milligrams
¹ Protein concentration is determined in a modified dye binding assay using Bovine Gamma Globulin (BGG) as a standard		

